

Product Datasheet

Canine TIMP1 ELISA Kit (Ready to Use) (orb549402)

Description	Canine Tissue Inhibitors Of Metalloproteinase 1 (TIMP1) Sandwich ELISA KitReady to Use: Pre-diluted ...
Reactivity	Canine
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P81546
Sensitivity	0.066ng/mL
Entrez	403816
Expiration Date	Please enquire.