

Product Datasheet

Canine TPS ELISA Kit (Ready to Use) (orb549407)

Description	Canine Tryptase (TPS) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Canine
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P15944
Sensitivity	13.3pg/mL
Entrez	100049001
Expiration Date	Please enquire.