

Product Datasheet

Gallus ADP ELISA Kit (Ready to Use) (orb549411)

Description	Chicken Adiponectin (ADP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni...
Reactivity	Gallus
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6QWE7
Sensitivity	0.292ng/mL
Entrez	404536
Expiration Date	Please enquire.