

Product Datasheet

Gallus CNTF ELISA Kit (Ready to Use) (orb549420)

Description	Chicken Ciliary Neurotrophic Factor (CNTF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection rea...
Reactivity	Gallus
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q02011
Sensitivity	29.4pg/mL
Entrez	396488
Expiration Date	Please enquire.