

Product Datasheet

Gallus C3 ELISA Kit (Ready to Use) (orb549422)

Description	Chicken Complement Component 3 (C3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents f...
Reactivity	Gallus
Range	31.2-2000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q90633
Sensitivity	11.2ng/mL
Entrez	396370
Expiration Date	Please enquire.