

Product Datasheet

Gallus FGF1 ELISA Kit (Ready to Use) (orb549425)

Description	Chicken Fibroblast Growth Factor 1, Acidic (FGF1) Sandwich ELISA Kit Ready to Use: Pre-diluted detect...
Reactivity	Gallus
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P19596
Sensitivity	21.54pg/mL
Entrez	396094
Expiration Date	Please enquire.