

Product Datasheet

Gallus GHRL ELISA Kit (Ready to Use) (orb549429)

Description	Chicken Ghrelin (GHRL) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient...
Reactivity	Gallus
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q7T2V1
Sensitivity	0.061ng/mL
Entrez	408185
Expiration Date	Please enquire.