

## Product Datasheet

### Gallus GC ELISA Kit (Ready to Use) (orb549430)

<b>Description</b>	Chicken Glucagon (GC) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
<b>Reactivity</b>	Gallus
<b>Range</b>	23.4-1500pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>P68259</b>
<b>Sensitivity</b>	7.25pg/mL
<b>Entrez</b>	<b>396196</b>
<b>Expiration Date</b>	Please enquire.