

## Product Datasheet

# Gallus HSPb2 ELISA Kit (Ready to Use) (orb549433)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Chicken Heat Shock Protein Beta 2 (HSPb2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent... |
| <b>Reactivity</b>        | Gallus   |
| <b>Range</b>             | 15.6-1000pg/mL   |
| <b>Form/Appearance</b>   | Ready to Use   |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt  |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.        |
| <b>Sample Types</b>      | tissue homogenates and other biological fluids.  |
| <b>Uniprot ID</b>        | <b>R4GG05</b>  |
| <b>Sensitivity</b>       | 5.7pg/mL   |
| <b>Entrez</b>            | <b>414136</b>  |
| <b>Expiration Date</b>   | Please enquire.  |