

## Product Datasheet

# Gallus IFN gamma ELISA Kit (Ready to Use) (orb549438)

<b>Description</b>	Chicken Interferon Gamma (IFN gamma) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents ...
<b>Reactivity</b>	Gallus
<b>Range</b>	15.6-1000pg/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>P49708</b>
<b>Sensitivity</b>	5.8pg/mL
<b>Entrez</b>	<b>396054</b>
<b>Expiration Date</b>	Please enquire.