

Product Datasheet

Gallus IL10 ELISA Kit (Ready to Use) (orb549440)

Description	Chicken Interleukin 10 (IL10) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Gallus
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q6A2H4
Sensitivity	32pg/mL
Entrez	428264
Expiration Date	Please enquire.