

Product Datasheet

Gallus IL12A ELISA Kit (Ready to Use) (orb549441)

Description	Chicken Interleukin 12A (IL12A) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Gallus
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q6X0L0
Sensitivity	7.2pg/mL
Entrez	407090
Expiration Date	Please enquire.