

Product Datasheet

Gallus MIP1b ELISA Kit (Ready to Use) (orb549451)

Description	Chicken Macrophage Inflammatory Protein 1 Beta (MIP1b) Sandwich ELISA Kit Ready to Use: Pre-diluted d...
Reactivity	Gallus
Range	62.5-4000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q90826
Sensitivity	25.8pg/mL
Entrez	395551
Expiration Date	Please enquire.