

Product Datasheet

Gallus MMP13 ELISA Kit (Ready to Use) (orb549453)

Description	Chicken Matrix Metalloproteinase 13 (MMP13) Sandwich ELISA Kit Ready to Use: Pre-diluted detection re...
Reactivity	Gallus
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O93363
Sensitivity	0.63ng/mL
Entrez	395683
Expiration Date	Please enquire.