

Product Datasheet

Gallus PGC ELISA Kit (Ready to Use) (orb549458)

Description	Chicken Pepsinogen C (PGC) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Gallus
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9PWK1
Sensitivity	0.198ng/mL
Entrez	395690
Expiration Date	Please enquire.