

Product Datasheet

Equine GH ELISA Kit (Ready to Use) (orb549473)

Description	Equine Growth Hormone (GH) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for conven...
Reactivity	Equine
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01245
Sensitivity	0.085ng/mL
Entrez	100034180
Expiration Date	Please enquire.