

Product Datasheet

Other ADP ELISA Kit (Ready to Use) (orb549521)

Description	General Adenosine Diphosphate (ADP) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents f...
Reactivity	All
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Sensitivity	0.101ng/mL
Expiration Date	Please enquire.