

Product Datasheet

Other ADMA ELISA Kit (Ready to Use) (orb549531)

Description	General Asymmetrical Dimethylarginine (ADMA) Sandwich ELISA KitReady to Use: Pre-diluted detection r...
Reactivity	All
Range	15.6-1000ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	4.63ng/mL
Expiration Date	Please enquire.