

Product Datasheet

Other E1 ELISA Kit (Ready to Use) (orb549553)

Description	General Estrone (E1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient a...
Reactivity	All
Range	125-8000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	62.4pg/mL
Expiration Date	Please enquire.