

Product Datasheet

Other GSH ELISA Kit (Ready to Use) (orb549561)

| | |
|--------------------------|--|
| Description | General Glutathione (GSH) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conveni... |
| Reactivity | All |
| Range | 1.56-100µg/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Sensitivity | 0.059µg/mL |
| Expiration Date | Please enquire. |