

Product Datasheet

Other LTB4 ELISA Kit (Ready to Use) (orb549569)

Description	General Leukotriene B4 (LTB4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	All
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.28ng/mL
Expiration Date	Please enquire.