

Product Datasheet

Other MT ELISA Kit (Ready to Use) (orb549572)

Description	General Melatonin (MT) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient...
Reactivity	All
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	6.5pg/mL
Expiration Date	Please enquire.