

Product Datasheet

Other NAADP ELISA Kit (Ready to Use) (orb549573)

| | |
|--------------------------|---|
| Description | General Nicotinic Acid Adenine Dinucleotide Phosphate (NAADP) Sandwich ELISA KitReady to Use: Pre-di... |
| Reactivity | All |
| Range | 0.156-10ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Sensitivity | 0.055ng/mL |
| Expiration Date | Please enquire. |