

Product Datasheet

Other NE ELISA Kit (Ready to Use) (orb549575)

Description	General Noradrenaline (NE) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	All
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	27.2pg/mL
Expiration Date	Please enquire.