

Product Datasheet

Other PAF ELISA Kit (Ready to Use) (orb549579)

Description	General Platelet Activating Factor (PAF) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reage...
Reactivity	All
Range	0.468-30ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.224ng/mL
Expiration Date	Please enquire.