

Product Datasheet

Other SCT ELISA Kit (Ready to Use) (orb549587)

Description	General Salmon Calcitonin (SCT) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for c...
Reactivity	All
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability biological agents.
Sample Types	biological agents.
Sensitivity	0.085ng/mL
Expiration Date	Please enquire.