

Product Datasheet

Other T ELISA Kit (Ready to Use) (orb549588)

Description	General Testosterone (T) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenience...
Reactivity	All
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, urine or other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Sensitivity	0.063ng/mL
Expiration Date	Please enquire.