

Product Datasheet

Guinea pig ALB ELISA Kit (Ready to Use) (orb549619)

Description	Guinea pig Albumin (ALB) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Guinea pig
Range	1.56-100µg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6WDN9
Sensitivity	0.48µg/mL
Entrez	100379232
Expiration Date	Please enquire.