

Product Datasheet

Guinea pig IFN gamma ELISA Kit (Ready to Use) (orb549622)

Description	Guinea pig Interferon Gamma (IFN gamma) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8CGS0
Sensitivity	6.3pg/mL
Entrez	100379568
Expiration Date	Please enquire.