

Product Datasheet

Human HTR4 ELISA Kit (Ready to Use) (orb549647)

Description	Human 5-Hydroxytryptamine Receptor 4 (HTR4) Sandwich ELISA Kit Ready to Use: Pre-diluted detection re...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q13639
Sensitivity	0.063ng/mL
Entrez	3360
Expiration Date	Please enquire.