

## Product Datasheet

### Human NT5E ELISA Kit (Ready to Use) (orb549650)

<b>Description</b>	Human 5'-Nucleotidase, Ecto (NT5E) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents fo...
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100ng/mL
<b>Form/Appearance</b>	Ready to Use
<b>Storage</b>	All reagents should be stored at -20°C upon receipt
<b>Note</b>	For research use only
<b>Application notes</b>	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>Q53Z63</b>
<b>Sensitivity</b>	0.66ng/mL
<b>Entrez</b>	<b>4907</b>
<b>Expiration Date</b>	Please enquire.