

Product Datasheet

Human ACTN3 ELISA Kit (Ready to Use) (orb549692)

Description	Human Actinin Alpha 3 (ACTN3) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for con...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q08043
Sensitivity	0.115ng/mL
Entrez	89
Expiration Date	Please enquire.