

Product Datasheet

Human ADCYAP1 ELISA Kit (Ready to Use) (orb549719)

Description	Human Adenylate Cyclase Activating Polypeptide 1, Pituitary (ADCYAP1) Sandwich ELISA KitReady to Use...
Reactivity	Human
Range	7.81-500pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P18509
Sensitivity	2.3pg/mL
Entrez	116
Expiration Date	Please enquire.