

Product Datasheet

Human ADIPOR2 ELISA Kit (Ready to Use) (orb549726)

Description	Human Adiponectin Receptor 2 (ADIPOR2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q86V24
Sensitivity	0.24ng/mL
Entrez	79602
Expiration Date	Please enquire.