

Product Datasheet

Human AFM ELISA Kit (Ready to Use) (orb549742)

Description	Human Afamin (AFM) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient and...
Reactivity	Human
Range	12.5-800ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P43652
Sensitivity	5.92ng/mL
Entrez	173
Expiration Date	Please enquire.