

Product Datasheet

Human a2PI ELISA Kit (Ready to Use) (orb549780)

Description	Human Alpha-2-Plasmin Inhibitor (a2PI) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	12.5-800ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P08697
Sensitivity	4.3ng/mL
Entrez	5345
Expiration Date	Please enquire.