

Product Datasheet

Human aEP ELISA Kit (Ready to Use) (orb549781)

Description	Human Alpha-Endorphin (aEP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell culture supernates or other biological fluid.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Sensitivity	4.45pg/mL
Expiration Date	Please enquire.