

Product Datasheet

Human APLP1 ELISA Kit (Ready to Use) (orb549807)

Description	Human Amyloid Beta Precursor Like Protein 1 (APLP1) Sandwich ELISA KitReady to Use: Pre-diluted dete...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P51693
Sensitivity	0.118ng/mL
Entrez	333
Expiration Date	Please enquire.