

Product Datasheet

Human ANXA10 ELISA Kit (Ready to Use) (orb549841)

Description	Human Annexin A10 (ANXA10) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Human
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9UJ72
Sensitivity	5.9pg/mL
Entrez	11199
Expiration Date	Please enquire.