

Product Datasheet

Human AP13 ELISA Kit (Ready to Use) (orb549856)

Description	Human Apelin 13 (AP13) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient...
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	37.2pg/mL
Expiration Date	Please enquire.