

Product Datasheet

Human AP36 ELISA Kit (Ready to Use) (orb549858)

Description	Human Apelin 36 (AP36) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient...
Reactivity	Human
Range	46.8-3000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	12.58pg/mL
Expiration Date	Please enquire.