

## Product Datasheet

# Human APEX1 ELISA Kit (Ready to Use) (orb549861)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Human APEX Nuclease 1 (APEX1)<br>Sandwich ELISA Kit<br>Ready to Use: Pre-diluted detection reagents for con...     |
| <b>Reactivity</b>        | Human  |
| <b>Range</b>             | 0.312-20ng/mL  |
| <b>Form/Appearance</b>   | Ready to Use   |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt  |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stability tissue homogenates, cell lysates or other biological fluids. |
| <b>Sample Types</b>      | tissue homogenates, cell lysates and other biological fluids.  |
| <b>Sensitivity</b>       | 0.132ng/mL   |
| <b>Expiration Date</b>   | Please enquire.  |