

Product Datasheet

Human AQP2 ELISA Kit (Ready to Use) (orb549883)

Description	Human Aquaporin 2, Collecting Duct (AQP2) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability tissue homogenates, urine or other biological fluids.
Sample Types	tissue homogenates, urine and other biological fluids.
Uniprot ID	P41181
Sensitivity	0.122ng/mL
Entrez	359
Expiration Date	Please enquire.