

Product Datasheet

Human ATM ELISA Kit (Ready to Use) (orb549923)

Description	Human Ataxia Telangiectasia Mutated (ATM) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q13315
Sensitivity	0.047ng/mL
Entrez	472
Expiration Date	Please enquire.