

Product Datasheet

Human ABCC8 ELISA Kit (Ready to Use) (orb549937)

Description	Human ATP Binding Cassette Transporter C8 (ABCC8) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q09428
Sensitivity	0.128ng/mL
Entrez	6833
Expiration Date	Please enquire.