

Product Datasheet

Human BAG3 ELISA Kit (Ready to Use) (orb549985)

Description	Human Bcl2 Associated Athanogene 3 (BAG3) Sandwich ELISA KitReady to Use: Pre-diluted detection reag...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	095817
Sensitivity	0.069ng/mL
Entrez	9531
Expiration Date	Please enquire.