

Product Datasheet

Human bACE1 ELISA Kit (Ready to Use) (orb549995)

| | |
|--------------------------|--|
| Description | Human Beta-Site APP Cleaving Enzyme 1 (bACE1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection ... |
| Reactivity | Human |
| Range | 0.25-16ng/mL |
| Form/Appearance | Ready to Use |
| Storage | All reagents should be stored at -20°C upon receipt |
| Note | For research use only |
| Application notes | Pre-diluted detection reagents and enhanced stability serum, plasma, cerebrospinal fluid or other biological fluids. |
| Sample Types | serum, plasma, cerebrospinal fluid and other biological fluids. |
| Uniprot ID | P56817 |
| Sensitivity | 0.11ng/mL |
| Entrez | 23621 |
| Expiration Date | Please enquire. |