

Product Datasheet

Human bTC ELISA Kit (Ready to Use) (orb549997)

Description	Human Betacellulin (bTC) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35070
Sensitivity	0.064ng/mL
Entrez	685
Expiration Date	Please enquire.