


Product Datasheet

Human Angiostatin Kringle 1-3/PLG ELISA Kit (orb55000)

Description	Human Angiostatin Kringle 1-3/PLG ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V98-S355 (K1-3)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1000, concentration in serum and plasma is around 1000 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
MW	90569 MW
Uniprot ID	P00747
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Angiostatin K1-3 PicoKine ELISA Ki...