

Product Datasheet

Human BST1 ELISA Kit (Ready to Use) (orb550006)

Description	Human Bone Marrow Stromal Cell Antigen 1 (BST1) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q10588
Sensitivity	33pg/mL
Entrez	683
Expiration Date	Please enquire.